


# CUT&RUN

---


The stronger version of ChIP-seq

Laureano Tomás-Daza  
SubGM3  
28/04/2021

# Cleavage Under Targets and Release Using Nuclease (CUT&RUN)




# CUT&RUN




- Add ConA-coated microscopic magnetic beads

# CUT&RUN




- Add ConA-coated microscopic magnetic beads
- They stick to the cell surface

# CUT&RUN



- Permeabilize plasmatic and nuclear membrane
- Add antibody, which is able to permeate inside the nucleus and bind


# CUT&RUN




- Permeabilize plasmatic and nuclear membrane
- Add antibody, which is able to permeate inside the nucleus and bind

# CUT&RUN


- Place on magnet and excess wash antibody



# CUT&RUN




- Remove from magnet
- We have intact permeabilized cells with the antibody bound to our protein of interest in native chromatin








# CUT&RUN




- Add pA/G-MNase fusion protein to the sample.
- p/AG binds (🏹) to the constant fraction of antibodies
- MNase (🌿) is an inespecific nuclease which cuts wherever it finds accessible

# CUT&RUN




- Add pA/G-MNase fusion protein to the sample.
- p/AG binds (  ) to the constant fraction of antibodies
- MNase (  ) is an inespecific nuclease which cuts wherever it finds accessible

# CUT&RUN




- Activation of Nuclease
- Nuclease cuts only in the regions where it has bound

# CUT&RUN




- Cut fragments are able to solubilize and permeate through the membrane; while genomic DNA does not


# CUT&RUN




- Cut fragments are able to solubilize and permeate through the membrane; while genomic DNA does not




# CUT&RUN



- Recover supernatant with our DNA fragments of interest, leaving behind the cells and the rest of the genomic DNA



# CUT&RUN



- Purify DNA

# CUT&RUN is ChIPseq on STEROIDS

The technique is performed on native chromatin, which ensures **preservation of the protein epitope**

CUT&RUN



No need for crosslinking  
or sonication

WAAAY less cells needed

Improved signal and less background

ChIPseq



Basically the same lame  
technique since 1985  
(Solomon and Varshavsky)



# CUT&RUN is ChIPseq on STEROIDS

Standard ChIPseq protocols  
require millions of cells;  
CUT&RUN developers claim it  
can work with works with as  
little as 100 cells for abundant  
histone marks or 1000 cells for  
TFs

CUT&RUN



No need for crosslinking  
or sonication

WAAAY less cells needed


Improved signal and less background

ChIPseq



Basically the same lame  
technique since 1985  
(Solomon and Varshavsky)

# CUT&RUN is ChIPseq on STEROIDS



*Adapted from Skene and Henikoff (2017)*

CUT&RUN



No need for crosslinking  
or sonication

WAAAY less cells needed


Improved signal and less background

ChIPseq




Basically the same lame  
technique since 1985  
(Solomon and Varshavsky)

# CUT&RUN Tools




Zhu, Q., Liu, N., Orkin, S.H. *et al.* (2019). CUT&RUNTools: a flexible pipeline for CUT&RUN processing and footprint analysis. *Genome Biol*

# Our (preliminary) CUT&RUN pipeline



# Our (preliminary) result




Have you thought "oh it seems really easy"?


**SPOILER:**

**It's not that easy**  
(yet)

# Critical Step 1: Alignment




# Critical Step 1: Alignment





## Critical Step 2: Peak Calling



# Critical Step 3: Normalization – Spike-in

## Normalization 1:

Treatment spike-in = 100,000 spike-in reads in 5,000,000 uniquely aligned reads = 2%

Untreated spike-in = 30,000 spike-in reads in 3,000,000 uniquely aligned reads = 1%

Ratio total/spike-in

Treatment normalization factor =  $1 / 2\%$  spike-in = **0.5**

Untreated normalization factor =  $1 / 1\%$  spike-in = **1.0**

## Normalization 2:

Treatment spike-in = 100,000 spike-in reads

Untreated spike-in = 30,000 spike-in reads


Absolute spike-in reads

Treatment normalization factor =  $30,000 / 100,000$  = **0.3**

Untreated normalization factor =  $30,000 / 30,000$  = **1.0**

# Critical Step 3: Normalization – Background

Suggested by: Samuel Collombet

1. Merge peak-calling from conditions
  2. Mask genome with 1.
  3. Coverage in 10Kb windows
  4. Compute scaling factor with DESeq2
- 
- Define background

# Take home message

- Both ChIPseq and CUT&RUN aim to answer the same biological question: where does a certain protein bind in the genome??

## HOWEVER

- CUT&RUN has lower cell input and sequencing requirements
- CUT&RUN gives lower background
- CUT&RUN maintains epitope integrity due to lack of crosslink and sonication

# Thanks for your attention

---

Any question?

Laureano Tomás-Daza  
SubGM3  
28/04/2021